



▶ Tri-anim®


ANESTHESIA PRODUCT CATALOG

2024

FEATURED PRODUCTS




p.4 **Curaplex® CuraView® LED Laryngoscopes**



p.10 **DART Intranasal Mucosal Atomization Device**



p.26 **Intersurgical® i-gel® Supraglottic Airway**



p.12 **RIVANNA® Accuro® Spinal Navigation System**



p.14 **Curaplex® PrO₂ mask**



p.31 **Curaplex® Accucuff™**



p.16 **Curaplex® EtCO₂ Sampling Cannula**



p.20 **Philips Nasal Alar SpO₂ Sensor by Curaplex®**

CONTENTS

- 4 Laryngoscope Blades & Handles
- 9 IV Products
- 9 Safety/Protection
- 10 Medication Delivery
- 12 Epidural Guidance
- 14 Capnography
- 20 Patient Monitoring
- 23 Anesthesia Masks
- 26 Supraglottic Airways
- 28 Endotracheal Tubes
- 32 Stylets
- 33 Introducers
- 34 Nasopharyngeal Airways
- 36 Oral Airways
- 39 Regional Anesthesia
- 40 Blenders & Hoses
- 41 CO₂ Absorbents
- 42 Adapters & Connectors
- 44 HME/Filters
- 47 Anesthesia Bags
- 48 Circuits



Tri-anim's Capital Rewards Program

Products marked with this icon are eligible for either our Rewards or Rewards Select program.

our brand

Curaplex® Brand

Curaplex® by Tri-anim is committed to helping you solve the unique challenges you face every day. Our anesthesia product solutions are inspired by and created for clinicians, to improve patient care while reducing overall healthcare costs.

Abbreviations

Fr - French | L - Length | W - Width | D - Diameter
ID - Inner Diameter | OD - Outer Diameter



CURAVIEW® LED FIBER OPTIC & LED ONE-PIECE LARYNGOSCOPES

The Curaplex® CuraView® LED system provides a completely sterile single-use laryngoscope system designed with industry-leading LED brightness for enhanced visibility for both everyday and difficult intubations. CuraView® combines the quality and strength of a reusable laryngoscope with the convenience and dependability of a disposable.

- Sterile**
Eliminate the risk of patient to patient cross-contamination and factory embedded contamination with a completely sterile single-use solution.
- Clear LED Illumination**
Regardless of the patient scenario and urgency, CuraView® LED provides clear LED lighting that you can count on for optimal visibility during intubation.
- No Hidden Cost**
Straightforward pricing eliminates the hidden cost associated with reusable laryngoscopes, such as replacement parts, cleaning supplies and labor.
- Quality**
The CuraView® LED is designed with a no-slip handle and low blade profile that inspires confidence during your everyday and emergent neonatal intubations.



TWO-PIECE OPTION
Curaplex® CuraView® LED Fiber Optic



ONE-PIECE OPTION
Curaplex® CuraView® LED One-Piece

our brand

Curaplex® CuraView® LED Fiber-Optic, Handles

Part #	Handle Type	Qty
301-F0-S	Standard	20/bx
301-F0-M	Mini	20/bx
301-F0-Y	Stubby	20/bx

our brand

Curaplex® CuraView® LED Fiber-Optic, Mac Blades

Part #	Size	Blade Type	Color Code	Qty
301-F0-MAC-1	1	Mac	Grey	20/bx
301-F0-MAC-2	2	Mac	Green	20/bx
301-F0-MAC-3	3	Mac	Orange	20/bx
301-F0-MAC-35	3.5	Mac	Purple	20/bx
301-F0-MAC-4	4	Mac	Red	20/bx

our brand

Curaplex® CuraView® LED Fiber-Optic, Miller Blades


Part #	Size	Blade Type	Color Code	Qty
NEW SIZE!				
301-F0-MIL-000	000	Miller	Dark Blue	20/bx
301-F0-MIL-00	00	Miller	Yellow	20/bx
301-F0-MIL-0	0	Miller	White	20/bx
301-F0-MIL-1	1	Miller	Grey	20/bx
301-F0-MIL-15	1.5	Miller	Pink	20/bx
301-F0-MIL-2	2	Miller	Green	20/bx
301-F0-MIL-3	3	Miller	Orange	20/bx
301-F0-MIL-4	4	Miller	Red	20/bx

our brand

Curaplex® CuraView® LED Fiber-Optic, Specialty Blades

Part #	Size	Blade Type	Color Code	Qty
301-F0-RS-0	0	Robert Shaw	White	20/bx
301-F0-PH-1	1	Phillips	Grey	20/bx
301-F0-RS-1	1	Robert Shaw	Grey	20/bx
301-F0-WH-1	1	Wis-Hipple	Grey	20/bx
301-F0-WH-15	1.5	Wis-Hipple	Pink	20/bx
301-F0-PH-2	2	Phillips	Green	20/bx
301-F0-WH-2	2	Wis-Hipple	Green	20/bx

LARYNGOSCOPE BLADES & HANDLES (SINGLE-USE)





our brand

Curaplex® CuraView® LED Fiber-Optic, Mac Combo Sets

Part #	Size	Blade Type	Handle Type	Color Code	Qty
301-FO-MAC-1-S	1	Mac	Standard		20/bx
301-FO-MAC-2-S	2	Mac	Standard		20/bx
301-FO-MAC-3-S	3	Mac	Standard		20/bx
301-FO-MAC-3.5-S	3.5	Mac	Standard		20/bx
301-FO-MAC-4-S	4	Mac	Standard		20/bx



our brand

Curaplex® CuraView® LED Fiber-Optic, Miller Combo Sets

Part #	Size	Blade Type	Handle Type	Color Code	Qty
301-FO-MIL-00-S	00	Miller	Standard		20/bx
301-FO-MIL-0-S	0	Miller	Standard		20/bx
301-FO-MIL-1.5-S	1.5	Miller	Standard		20/bx
301-FO-MIL-1-S	1	Miller	Standard		20/bx
301-FO-MIL-2-S	2	Miller	Standard		20/bx
301-FO-MIL-3-S	3	Miller	Standard		20/bx
301-FO-MIL-4-S	4	Miller	Standard		20/bx



our brand

Curaplex® CuraView® LED Fiber-Optic, Mini Handle Combo Sets

Part #	Size	Blade Type	Handle Type	Color Code	Qty
NEW SIZE! 301-FO-MIL-000-M	000	Miller	Mini		20/bx
301-FO-MIL-00-M	00	Miller	Mini		20/bx
301-FO-MIL-0-M	0	Miller	Mini		20/bx
301-FO-MAC-1-M	1	Mac	Mini		20/bx
301-FO-MIL-1-M	1	Miller	Mini		20/bx
301-FO-MIL-1.5-M	1.5	Miller	Mini		20/bx



our brand

Curaplex® CuraView® LED Fiber-Optic, Specialty Combo Sets

Part #	Size	Blade Type	Handle Type	Color Code	Qty
301-FO-PH1-S	1	Phillips	Standard		20/bx
301-FO-PH2-S	2	Phillips	Standard		20/bx

our brand

Curaplex® CuraView® LED, Mac Blades

Part #	Size	Blade Type	Handle Type	Color Code	Qty
301-LED-010S	1	Mac	Standard		20/bx
301-LED-020S	2	Mac	Standard		20/bx
301-LED-030S	3	Mac	Standard		20/bx
301-LED-035S	3.5	Mac	Standard		20/bx
301-LED-040S	4	Mac	Standard		20/bx
301-LED-010M	1	Mac	Mini		20/bx
301-LED-020M	2	Mac	Mini		20/bx

our brand

Curaplex® CuraView® LED, Miller Blades

Part #	Size	Blade Type	Handle Type	Color Code	Qty
301-LED-100S	00	Miller	Standard		20/bx
301-LED-000S	0	Miller	Standard		20/bx
301-LED-101S	1	Miller	Standard		20/bx
301-LED-1015S	1.5	Miller	Standard		20/bx
301-LED-102S	2	Miller	Standard		20/bx
301-LED-103S	3	Miller	Standard		20/bx
301-LED-104S	4	Miller	Standard		20/bx
NEW SIZE! 301-LED-1000M	000	Miller	Mini		20/bx
301-LED-100M	00	Miller	Mini		20/bx
301-LED-000M	0	Miller	Mini		20/bx
301-LED-101M	1	Miller	Mini		20/bx
301-LED-1015M	1.5	Miller	Mini		20/bx

our brand

Curaplex® CuraView® LED, Specialty Blades

Part #	Size	Blade Type	Handle Type	Color Code	Qty
301-LED-PH1-S	1	Phillips	Standard		20/bx
301-LED-WH1-S	1	Wis-Hipple	Standard		20/bx
301-LED-WH15-S	1.5	Wis-Hipple	Standard		20/bx
301-LED-PH2-S	2	Phillips	Standard		20/bx
301-LED-WH2-S	2	Wis-Hipple	Standard		20/bx
301-LED-PH1-M	1	Phillips	Mini		20/bx
301-LED-WH1-M	1	Wis-Hipple	Mini		20/bx
301-LED-WH15-M	1.5	Wis-Hipple	Mini		20/bx
301-LED-PH2-M	2	Phillips	Mini		20/bx





our brand

Curaplex® Reusable Blades, FO Green System

Item #	Size	Handle Type	Style	Qty
301-B4011	1	Infant	Mac	1/ea
301-B4012	2	Child	Mac	1/ea
301-B4013	3	Medium Adult	Mac	1/ea
301-B4013S	3	Strong Curve	Mac	1/ea
301-B4014	4	Large Adult	Mac	1/ea
301-B40200	00	Premature	Miller	1/ea
301-B4020	0	Neonate	Miller	1/ea
301-B4021	1	Infant	Miller	1/ea
301-B4022	2	Child	Miller	1/ea
301-B4023	3	Medium Adult	Miller	1/ea
301-B4024	4	Large Adult	Miller	1/ea



our brand

Curaplex® Reusable Handles, FO Green System with Halogen

Item #	Handle Type	Qty
301-B701	Standard	1/ea
301-B702	Mini	1/ea

our brand

Curaplex® Reusable Handles, FO Green System with Premium LED*

Item #	Handle Type	Qty
301-B850	Mini	1/ea
301-B851	Stubby	1/ea
301-B852	Standard	1/ea
301-B853	Premium Mini	1/ea
301-B854	Premium Standard	1/ea

*Premium LED handles have a lifetime warranty

our brand

Curaplex® Pressure Infusers

Pressure pop-off activates at 325 mmHg.

Item #	Volume	Qty
301-MTM-310EA	1000 ml	25/ca
301-MTM305EA	500 ml	25/ca



AirLife® Pressure Infusion Bag

AirLife Pressure Infusion bags are made of durable vinyl and feature a color-coded pressure gauge for easy classification and monitoring. The see-through mesh material allows for easy visibility and increased safety.

Item #	Size	Qty
87-VC430EA	3000 ml	10/bx
87-VC410EA	1000 ml	25/bx
87-VC405EA	500 ml	25/bx



our brand

Curaplex® Eye Protector

Curaplex® eye protectors are rigid, clear, plastic protective covers used to shield the patients’ eyes during surgical procedures. The eye protectors feature a non-allergenic, self-adhesive foam cushion to enable quick and accurate application on the patient. The thick foam padding increases patient comfort even when the patient is in the prone position. The one-piece backing has “courtesy tabs” to assist with removal. Latex-free, sterile and single-use. Made in the USA.

Item #	Size	Qty
301-AEG100EA	Adult	25/bx



AirLife® IGuard™ Eye Protector

- Protects patients eyes during general anesthesia against eye abrasion, lacerations and scleral hemorrhage
- Clear plastic protective covers with non-allergenic self-adhesive foam
- Improved tab design means one-handed release of adhesive paper
- Thicker adhesive foam is more comfortable for patient
- Fast and accurate application, individually packaged
- Disposable, latex-free, sterile

Item #	Size	Qty
792-9-0210-00	Adult	25/ca
792-9-0210-20	Adult, Double Foam	25/ca
792-9-0210-21	Pediatric	25/ca





DART INTRANASAL MUCOSAL ATOMIZATION DEVICE

Deliver medications quickly and painlessly with the DART. Atomized medication delivery results in a therapeutic onset of action similar to intravenous delivery - with no pain or delay associated with starting an IV.


Features:

- Contain soft conical-shaped foam tip to facilitate comfortable application and forms a good seal within the nostril
- Measure and deliver with accuracy
- No IV line needed

our brand

Curaplex® DART Intranasal Mucosal Atomization Device

Item #	Description	Qty
301-2170-20300	Curaplex® DART without syringe	25/ca
301-2170-20150	Curaplex® DART with 1mL syringe	25/ca
301-2170-20100	Curaplex® DART with 3mL syringe	25/ca



DART REACH LARYNGO-TRACHEAL MUCOSAL ATOMIZATION DEVICE

Deliver medications quickly and effectively for difficult airways with DART Reach. Features a flexible stylet that gives control over the procedure and the ability to adjust the tip based on the patient’s anatomy.


Features:

- Flexible stylet that can be adjusted based on the patient’s anatomy
- ISO standard luer lock
- Available in adult and pediatric sizes

our brand

Curaplex® DART Reach Laryngo-Tracheal Musocal Atomization Device

Item #	Description	Qty
301-REACH720	Curaplex® DART Reach (4 inches) Device without 5mL Syringe	25/ca
301-REACH700	Curaplex® DART Reach (8 inches) Device without 5mL Syringe	25/ca
301-REACH600	Curaplex® DART Reach (8 inches) Device with 5mL Syringe	25/ca



Aerogen® Solo

Take the guess work out of aerosol medication delivery. Aerogen Solo offers a hands-free, high performance aerosol medication delivery system which provides optimal patient care.¹

Item #	Description	Qty
06-AG-AS7500-US	Aerogen Ultra Starter Kit	1/ea
06-AG-AS4300-US	OR Convenience Kit	1/ea
06-AG-PX1000-US	Aerogen Pro-X Controller	1/ea
06-AG-AS7000-US	Aerogen Ultra	1/ea
203-R100P90	Maxtec Stand and Basket	1/ea
203-R100P68	Maxtec Aerogen Locking Bracket	1/ea



PROTECTION FOR PATIENTS & CAREGIVERS

STERILE AND SINGLE-USE



Curaplex® CuraView® LED Fiber-Optic Laryngoscope

HIGH EFFICIENCY SAFETY



Intersurgical® Filta-Therm™ Heat and Moisture Exchanging Filters

CLOSED SYSTEM THAT WON'T BREAK THE CIRCUIT



Aerogen® Solo Convenience Kits

RIVANNA®

CERTAINTY CAN
BE EFFORTLESS

RIVANNA® Accuro® Spinal Navigation System

The Accuro® spinal navigation system delivers automated, real-time image guidance with significantly enhanced bone-to-tissue contrast to ensure improved localization of the desired intervertebral space for first-attempt success during spinal needle guidance procedures.



Hit the spot without a steep learning curve.

- 1. IDENTIFY** epidural location automatically with success rates exceeding 94%* using SpineNav3D™ AI-Based Spine Recognition.
- 2. MARK** needle placement or perform real-time needle guidance with the Accuro® Locator needle guide.
- 3. FIND** the ideal insertion point with Midline (red dashed line) and Cross Hair indicators.
- 4. SEE** more than 5 to 10X* enhancement of bone-to-tissue contrast with Multi-Frequency BoneEnhance® Image Reconstruction.

RIVANNA® Accuro® Epidural and Spinal Needle Guidance Instrument

Accuro® guides you to the desired intervertebral space with easy-to-interpret graphical indicators and enhanced bony landmark visualization.

Item #	Description	Service Plan	Qty
740-901-00002-3YEAR	Accuro® Epidural and Spinal Needle Guidance Instrument	3-Year	1/ea
740-901-00002	Accuro® Epidural and Spinal Needle Guidance Instrument	—	1/ea
740-901-00002-1YEAR	Accuro® Epidural and Spinal Needle Guidance Instrument	1-Year	1/ea
740-901-00002-2P-TA	Accuro® Epidural and Spinal Instrument, 2 Device Package	3-Year	1/ea
740-901-00002-CP-1YEAR	Certified Pre-Used Accuro® Epidural and Spinal Needle Guidance Instrument	1-Year	1/ea
740-901-00002-HSP-TA	Accuro® Head Start Package including 3 Accuro® Instruments, 3 Accuro Trolleys	3-Year	1/ea




740-901-00002

Accuro® Locator Needle Guide

The Locator is a patented single-use needle guide used with Accuro® to mark the needle insertion site accurately along the skin surface or perform real-time needle guidance.

Item #	Description	Qty
740-4742	Locator Needle Guide	24/bx



740-4742

RIVANNA® Sterile Cover and Locator Kit

The Locator and sterile cover contribute to a reduced risk of infection and patient motion during pre-procedural imaging, and allow for a re-scan after initial needle advancement in challenging patients if needed.

Item #	Description	Qty
740-4741	Kit containing Soft Pullup™ instrument cover, Locator needle guide, sterile gel, and 2 bands	24/bx



740-4741

Accuro® Trolley

The Accuro® sits upright in the top basket—under lock and key—so you can securely park your trolley on the hospital floor, ensuring your Accuro is safe and easily accessible. The 10,000 mAh power bank integrated under the Trolley lid makes it easy to keep the Accuro® charged-up and ready to navigate.

Item #	Description	Qty
740-4300	Trolley	1/ea




740-4300



PRO₂ MASK

Prevent hypoxemia with the Curaplex® PrO₂ Mask — Pro-actively support oxygenation before, during, and after procedures.

- Supports adequate oxygen delivery throughout procedure
- Dual oral and nasal access ports
- Provides higher FiO₂ than standard EtCO₂ cannula
- Mask size and shape enhances patient comfort
- Monitors continuous waveform capnography
- Single device for medium and high oxygen requirements



301-PrO2LT



our brand

Curaplex® PrO₂ Mask

High concentration (10-15 L/MIN ~ 86%-97% O₂) delivers high FiO₂ to your endoscopy patient

Item #	Description	Mask & O ₂ line	Oxygen Reservoir Bag	Medium Concentration Elbow	Microstream Sampling Line	M2M Sampling Line	M2F Sampling Line	Qty
301-PrO2LT	Adult O ₂ Mask - EtCO ₂ - M2M	✓	✓	✓		✓		30/bx
301-PrO2LT-MSA	Adult O ₂ Mask - EtCO ₂ - MSA	✓	✓	✓	✓			30/bx
301-PrO2LT-FM	Adult O ₂ Mask - EtCO ₂ - M2F	✓	✓	✓			✓	30/bx
301-PrO2LT-0	Adult O ₂ Mask - EtCO ₂	✓	✓	✓				30/bx



301-PrO2LTEZ-MM



301-PrO2LTEZ

Also available in:



MSA Sampling Line
301-PrO2LT-MSA



M2F Sampling Line
301-PrO2LT-FM

Also available in:



MSA Sampling Line
301-PrO2LTEZ-MSA



M2F Sampling Line
301-PrO2LTEZ-MF

Medium concentration (8-12 L/MIN ~ 81%- 95% O₂) supports intermediate levels of supplemental O₂ delivery

Item #	Description	Mask & O ₂ line	Oxygen Reservoir Bag	Medium Concentration Elbow	Microstream Sampling Line	M2M Sampling Line	M2F Sampling Line	Qty
301-PrO2LTEZ-MM	Adult O ₂ Mask - EtCO ₂ - M2M	✓		✓		✓		30/bx
301-PrO2LTEZ-MSA	Adult O ₂ Mask - EtCO ₂ - MSA	✓		✓	✓			30/bx
301-PrO2LTEZ	Adult O ₂ Mask - EtCO ₂	✓		✓				30/bx
301-PrO2LTEZ-MF	Adult O ₂ Mask - EtCO ₂ - M2F	✓		✓			✓	30/bx



ETCO₂ SAMPLING CANNULAS

EtCO₂/O₂ sampling cannulas collect exhaled CO₂ while simultaneously delivering oxygen. These EtCO₂/O₂ cannulas also feature a reflective style connector, compatible with Microstream™, Capnostream™, and Oridion™ capnography monitors.



In-line hydrophobic filter for improved moisture management

The tubing and ear piece are designed to help prevent kinking

Curved nares made of ultra-soft material

Fits-All O₂ connector comes standard on all sizes and lengths

Reflective style connector with ergonomic grip

Soft-Ears™ material helps ensure superior patient comfort, even during prolonged use


Also available in oral/nasal configuration



our brand

Standard EtCO₂ Cannulas

Standard use dual cannulas available in Nasal and Oral/Nasal configuration.



Orange reflective connector indicates the standard cannulas



Nasal Cannula




Oral/Nasal Cannula



our brand

Extended Use EtCO₂ Cannulas

Extended use cannulas are fitted with a section of Nafion tubing that allows the cannula to function longer and in higher-humidity conditions.



Yellow reflective connector indicates the extended line with Nafion tubing



Nasal Cannula



Oral/Nasal Cannula




our brand

Sampling Line with Airway Adapter (T-piece)

The Sampling line with Airway Adapter (T-Piece) can be used on both intubated and non-intubated patients. Available in both standard and extended use for Adult or Pediatric patients with ET tube size > 4.5 mm.



Standard Use



Extended Use

Item #	Description	Configuration			Size		Connector/Length of Use		Length of Sample Line			Qty
		Nasal Cannula	Oral/Nasal Cannula	Airway Adapter	Adult	Pediatric	Standard (Orange)	Extended use (Yellow)	7'	10'	14'	
301-5107F-SE	Adult EtCO ₂ /O ₂ Nasal Cannula	✓			✓		●		✓			25/bx
301-5110F-SE	Adult EtCO ₂ /O ₂ Nasal Cannula	✓			✓		●			✓		25/bx
301-5114F-SE	Adult EtCO ₂ /O ₂ Nasal Cannula	✓			✓		●				✓	25/bx
301-5107F-SE-E	Adult Extended EtCO ₂ /O ₂ Nasal Cannula	✓			✓			●	✓			25/bx
301-5114F-SE-E	Adult Extended EtCO ₂ /O ₂ Nasal Cannula	✓			✓			●			✓	25/bx
301-5207F-SE	Pediatric EtCO ₂ /O ₂ Nasal Cannula	✓				✓	●		✓			25/bx
301-5214F-SE	Pediatric EtCO ₂ /O ₂ Nasal Cannula	✓				✓	●				✓	25/bx
301-5207F-SE-E	Pediatric Extended EtCO ₂ /O ₂ Nasal Cannula	✓				✓		●	✓			25/bx
301-5214F-SE-E	Pediatric Extended EtCO ₂ /O ₂ Nasal Cannula	✓				✓		●			✓	25/bx
301-5707F-SE	Adult EtCO ₂ /O ₂ Oral/Nasal Cannula		✓		✓		●		✓			25/bx
301-5714F-SE	Adult EtCO ₂ /O ₂ Oral/Nasal Cannula		✓		✓		●				✓	25/bx
301-5707F-SE-E	Adult Extended EtCO ₂ /O ₂ Oral/Nasal Cannula		✓		✓			●	✓			25/bx
301-5714F-SE-E	Adult Extended EtCO ₂ /O ₂ Oral/Nasal Cannula		✓		✓			●			✓	25/bx
301-5107ET	Adult/Pediatric Sampling Line w/ Airway Adapter Intubated			✓	✓	✓	●		✓			25/bx
301-5107ET-E	Adult/Pediatric Sampling Line w/ Airway Adapter Intubated			✓	✓	✓		●	✓			25/bx
301-5114ET	Adult/Pediatric 14FT Sampling Line w/ Airway Adapter Intubated			✓	✓	✓	●				✓	25/bx

Microstream™, Capnostream™, and Oridion™ are trademarks of Medtronic, related to Product marked by Medtronic.




AirLife® EtCO₂ Scoop Cannulas

AirLife's EtCO₂ cannula simultaneously delivers oxygen while capturing CO₂. A series of strategically positioned holes deliver O₂ while CO₂ is effectively captured through the nares and scoop. This design helps lessen attenuation of the end tidal sample by oxygen dilution, particularly for mouth-breathing patients. The EtCO₂ cannula is non-sterile, latex-free and DEHP-free.

- Designed to prevent mixing of fresh oxygen with CO₂
- Soft Ears™ material ensures superior patient comfort, even during prolonged use
- 3-Channel tubing prevents kinking
- Includes male/female adapter with luer lock connector

Item #	Description	Length	Qty
87-4707	Adult EtCO ₂ Cannula	7' Tubing	25/ca
87-4710	Adult EtCO ₂ Cannula	10' Tubing	25/ca
87-4714	Adult EtCO ₂ Cannula	14' Tubing	25/ca



87-4707

Intersurgical Senti™ EcoLite™ EtCO₂ Mask

The Senti Mask is an adult, medium concentration EcoLite™ oxygen mask with a Luer lock connector enabling the connection of a monitoring line, the other end of which is connected to a capnometer. The Senti mask is used for measuring end-tidal CO₂ in expired gases while supplying oxygen to the patient.

Item #	Description	Qty
441-11416500	Adult, mask with CO ₂ monitoring line for Microstream® capnography with filter and tube, 2.1m	25/ca
441-1141015	Adult, mask with tube, 2.1m	30/ca
441-1142015	Adult, mask with CO ₂ monitoring line and tube, 2.1m	30/ca
441-11426100	Adult, mask with CO ₂ monitoring line with Universal Connector tube, 1.8m	30/ca



441-1142015



AirLife® EtCO₂ Divided Cannulas

AirLife's design provides true EtCO₂ sampling and proven O₂ delivery without limitations at higher O₂ rates. The divided design separates gases and Salter Eyes ensure an uninterrupted reading. OralTrac technology provides nasal sampling and an adjustable oral sampling trunk.

Item #	Description	Length	Luer Type	Qty
77-4796F-25	OralTrac, Adult, 2" CO ₂ Pigtail	7' Tubing	Female	25/ca
77-4797-25	OralTrac, Adult, Adapter	7' Tubing	Male	25/ca
77-4797F-25	OralTrac, Adult, Adapter	7' Tubing	Female	25/ca



77-4797F-25

AirLife® CapnoVue Mask

CapnoVue capnography masks improve capnography with oral nasal sampling for mouth breathers and upper GI cases and higher oxygen delivery for patients where more aggressive oxygenation is needed. Luer ports located at both expiratory ports separate CO₂ from O₂ for clean signals.

Item #	Description	Length	Luer Type	Qty
77-4727-7-0-50	Adult EtCO ₂ Mask	7' Tubing	Male	50/ca



77-4727-7-0-50

NASAL ALAR SPO₂ SENSOR BY CURAPLEX®

The Nasal Alar SpO₂ Sensor is your "one-and-done" solution for fast, accurate and dependable SpO₂ readings.



Optimal Site

Nasal alar is the highest perfused site with easy access during surgery.¹



Faster

Earlier detection 30 seconds faster* (2,3)
*compared to finger sensor



Reliable

Provides a signal even when other sensors fail.¹



Accurate

Accurate SpO₂ monitoring even with poorly perfused patients.



The Nasal Alar is rich in vasculature, making it an optimal pulse oximetry monitoring site. This central site provides faster detection of SpO₂ changes than the finger and offers significantly less signal loss due to poor perfusion.

1. Schallom M, Prentice D, Sona C, Arroyo C, Mazuski J. Comparison of nasal and forehead oximetry accuracy and pressure injury in critically ill patients. *Heart & Lung*. 2018 Mar 1; 47(2): 93-9.
2. Meiker, RJ, et al. Earlier Detection of Desaturation NADIR from the Nasal Ala. Presented at the Annual Meeting of the Society for Technology in Anesthesia, January 2013.
3. Comparative Desaturation Study performed at Xhale Assurance, Inc. and University of Florida, 2012. Data on file at Xhale Assurance, Inc.



Philips Nasal Alar SpO₂ Sensor by Curaplex®

As your pulse oximetry solution across the hospital, this single-patient-use sensor maintains placement on the nasal ala comfortably, without adhesives. The sensor is easily removed and reapplied, eliminating the waste associated with trying multiple sensors to get a good signal.

- Seven days of continuous use
- Dependable monitoring
- Earlier detection of SpO₂ changes
- Easy access to monitoring site
- Adhesive free

Item #	Description	Qty
301-11214	Philips Nasal Alar SpO ₂ ™ Sensor by Curaplex®	24/ca
73-M1020B	Philips FAST SpO ₂ Plug-in Module	1/ea
73-M1943AL	Philips Compatible SpO ₂ Adapter Cable	1/ea
73-989803205381	Nasal Alar FAST SpO ₂ Sensor	20/ca

RELIABLE OXIMETRY FOR POORLY PERFUSED PATIENTS

73-M1020B



Philips FAST Oximetry Module

73-98980320538



Philips FAST Nasal Alar Sensor

Earn a Philips FAST Oximetry Module for critical care patients in your facility through existing and future Philips FAST Nasal Alar Sensor purchases.

Tri-anim understands the business of healthcare and offers value-oriented programs designed to improve patient care and maximize operational efficiencies. The Rewards Select program provides immediate use of the Philips FAST Oximetry Module in return for consistent annual purchases of the Philips FAST Nasal Alar Sensors.

PROGRAM BENEFITS:

- ✓ Extend capital budget and comply with federal "Safe Harbor" laws
- ✓ Implement leading-edge technology into your facility to support patient care
- ✓ Receive capital equipment up front when anticipated disposable purchases are made



For more information, contact your dedicated Account Manager or learn more at: tri-anim.com/thc/capital-rewards



203-R113P90



203-R107P20



203-R112P10



For a complete listing of oxygen sensors, please contact your Tri-anim Account Manager.



Maxtec® Oxygen Sensors

The comprehensive selection of Maxtec® oxygen sensors can fulfill virtually all of your anesthesia applications.

Item #	Description	Qty
203-R113P90	Oxygen Sensor, MAX-1	1/ea
203-R107P20	Oxygen Sensor, MAX-2	1/ea
203-R107P22	Oxygen Sensor, MAX-8	1/ea
203-R112P10	Oxygen Sensor, MAX-10	1/ea
203-R113P10	Oxygen Sensor, MAX-11	1/ea



Maxtec MaxO₂ ME Monitor

MaxO₂ ME is an easy-to-use, low maintenance oxygen monitor. New features include a back lit LCD, smart alarms, DC power port, protective over-mold, built in kickstand, built in dove tail and extended battery life.

- Large, Back-Lit LCD
- Smart Alarms - auto set alarms 3% +/- current reading
- 15% low to 99% low alarm capabilities
- Long Battery Life - Up to 5,000 hours on 4 AA batteries
- DC Power Port, for use with optional remote power supply
- External Max-550E O₂ Sensor
- Built in kickstand and dovetail

Item #	Description	Qty
203-R230P01	MaxO ₂ ME Monitor	1/ea



Maxtec MaxO₂ ME+p Monitor

Continuously monitor oxygen and pressure during bubble CPAP therapy with the MaxO2 ME+p. Device can be integrated into a blender delivery system.

- Moisture-controlling Pressure Monitoring Line to connect to bCPAP system - unique Nafion™ Tubing Technology controls moisture in line

Item #	Description	Qty
203-R230P06	MaxO2 ME+p Oxygen and Pressure Monitor	1/ea



87-VR2600



87-VR3600



87-VR1600EA

AirLife® Soft Seal™ Mask

Ventlab’s disposable, non-inflatable anesthesia masks are designed specifically for pediatrics and are available in a variety of sizes.

Item #	Size	Qty
87-VR1600EA	Neonatal	20/ca
87-VR2600	Infant	20/ca
87-VR3600	Child	20/ca

SunMed® Vent Mask™ (with Valve)

The disposable Vent Mask™ is a premium inflatable, soft-cushioned pliable mask that has an ultra thin wall and an air-filled cushion to provide a low-pressure soft seal and superior seating characteristics.

Item #	Size	Qty
87-VR1100EA	Small Infant	20/ca
87-VR1100-40	Small Infant	40/ca
87-VR2100EA	Infant	20/ca
87-VR2100-40	Infant	40/ca
87-VR3100EA	Small Child	20/ca
87-VR3100-40	Small Child	40/ca
87-VR4100EA	Child	20/ca
87-VR4100-40	Child	40/ca
87-VR5000-20	Adult, Wide	20/ca
87-VR5000-50	Adult, Wide	50/ca
87-VR5000-100	Adult, Wide	100/ca
87-VR5100EA	Adult, Regular	20/ca
87-VR5100-50	Adult, Regular	50/ca
87-VR5100-100	Adult, Regular	100/ca
87-VR6100EA	Adult, Large	20/ca
87-VR6100-50	Adult, Large	50/ca

87-VR6100EA

87-VR5000-20

87-VR5100EA



87-VR4100EA



87-VR3100EA



87-VR2100EA



87-VR1100EA



AirLife® Soft Aire™ Mask (with Valve)

This disposable, inflatable, soft-cushion mask offers the same superior comfort and seal as other Ventlab™ masks and also includes a valve.

Item #	Size	Qty
87-VR1700EA	Small Infant	20/ca
87-VR2700EA	Infant	20/ca
87-VR2700-40	Infant	40/ca
87-VR2700-100	Infant	100/ca
87-VR3700-20	Small Child	20/ca
87-VR4700-20	Child	20/ca
87-VR5700-20	Adult	20/ca
87-VR5700-50	Adult	50/ca
87-VR5700-100	Adult	100/ca
87-VR6700-20	Adult, Wide	20/ca
87-VR6700-50	Adult, Wide	50/ca

Intersurgical® Economy Anesthetic Face Masks (with Color-coded Hook Rings)

- Color-coded hook ring: for ease of identification, and integration with other Intersurgical airway products
- Clear flexible housing and cushion: provides great patient visibility
- Soft cushion: for patient comfort and an air tight fit
- Four sizes available: from small adult to extra large adult

Item #	Size	Qty
441-1514000	3, Small Adult	50/bx
441-1515000	4, Medium Adult	40/bx
441-1516000	5, Large Adult	30/bx
441-1517000	6, Extra Large Adult	25/bx



441-1517000



AirLife® Vent Mask II™

The disposable Vent Mask II™ is a non-inflatable soft cushion mask that offers the same comfort and superior seal as the original Vent Mask™.

Item #	Size	Qty
87-VR0040EA	X-Small Neonatal	20/ca
87-VR0040-40	X-Small Neonatal	40/ca
87-VR0400EA	Small Neonatal	20/ca
87-VR0400-40	Small Neonatal	40/ca
87-VR1400-20	Small Infant	20/ca
87-VR1400EA	Small Infant	40/ca
87-VR2400-20	Infant	20/ca
87-VR2400EA	Infant	40/ca
87-VR2400-100	Infant	100/ca
87-VR3400-20	Small Child	20/ca
87-VR3400-40	Small Child	40/ca
87-VR4400EA	Child	20/ca
87-VR4400-40	Child	40/ca
87-VR5400EA	Adult	20/ca
87-VR6400EA	Large Adult	20/ca

Intersurgical® ClearLite™ Masks

ClearLite offers a range of anesthetic face masks which include seven sizes to meet the clinical needs of most patient types and sizes. Due to the enhanced anatomical design, the ClearLite mask offers superior seal performance with minimal pressure, reducing the risk of unnecessary discomfort or potential trauma. Crystal clear body provides excellent visibility for confirmation that the patient is well perfused. ClearLite masks are non-PVC and phthalate free which significantly minimizes the impact on the environment.

Item #	Description	Size	Qty
441-7290000	Blue Seal, No Hook Ring, 150D	0, Neonate	40/bx
441-7291000	Grey Seal, No Hook Ring, 150D	1, Infant	30/bx
441-7292000	White Seal, 22ID	2, Pediatric	25/bx
441-7293000	Yellow Seal, 22ID	3, Small Adult	35/bx
441-7294000	Green Seal, 22ID	4, Adult	35/bx
441-7295000	Orange Seal, 22ID	5, Large Adult	35/bx
441-7296000	Red Seal, 22ID	6, Extra Large Adult	20/bx



441-7296000

INTERSURGICAL®
I-GEL®
SUPRAGLOTTIC AIRWAY

The i-gel accurately and naturally positions itself over the laryngeal anatomy, providing a reliable perilaryngeal seal without the need for an inflating cuff.



441-8203030


- ✓

NON-INFLATING CUFF
 Made from a unique, soft, gel-like material allowing ease of insertion and reduced trauma
- ✓

INTEGRATED BITE BLOCK
 Reduces the possibility of airway channel occlusion
- ✓

GASTRIC CHANNEL
 Allows for the passing of a nasogastric tube to empty the stomach contents and facilitates venting (except size 1)
- ✓

CONDUIT FOR INTUBATION
 Under fiberoptic guidance in a known or unexpectedly difficult intubation (adult sizes: 3, 4, 5)



Intersurgical® i-gel®


The innovative second generation of supraglottic airway devices from Intersurgical. The first major development since the laryngeal mask airway, i-gel has changed the face of airway management and is now widely used in anesthesia and resuscitation across the globe.

Item #	Size	Qty
441-8201030	1, Neonatal	10/ca
441-8215030	1.5, Infant	10/ca
441-8202030	2, Small Pediatric	10/ca
441-8225030	2.5, Large Pediatric	10/ca
441-8203030	3, Small Adult	25/ca
441-8204030	4, Medium Adult	25/ca
441-8205030	5, Large Adult	25/ca

Intersurgical® Solus™ Standard Laryngeal Mask

Intersurgical Solus™ Laryngeal Masks are available in a wide range of comprehensive, disposable options including single-use Standard, Satin and Flexible Wire-reinforced.

Item #	Size	Qty
441-8001	1	20/bx
441-8002	2	20/bx
441-8003	3	20/bx
441-8003-F	3, Flexible	20/bx
441-8004	4	20/bx
441-8004-F	4, Flexible	20/bx
441-8005	5	20/bx
441-8005-F	5, Flexible	20/bx



441-8004

441-8004-F

Intersurgical® Solus™ Satin Laryngeal Mask

Solus Satin is our range of single-use laryngeal mask airways featuring a softer airway tube providing increased flexibility.

Item #	Size	Qty
441-8003003	3, Small Adult	20/bx
441-8004003	4, Medium Adult	20/bx
441-8005003	5, Large Adult	20/bx



441-8004003

Intersurgical® Solus™ MRI-Safe Laryngeal Mask

A range of single-use, sterile, laryngeal mask airways fitted with specially designed non-ferrous valves, guaranteed not to interfere with the magnet in an MRI scanner. The plastic valve has been selected to ensure full reliability throughout the shelf life of each device.

Item #	Size	Qty
441-8001002	1, Neonate	20/bx
441-8015002	1.5, Infant	20/bx
441-8002002	2, Small Pediatric	20/bx
441-8025002	2.5, Large Pediatric	20/bx
441-8003002	3, Small Adult	20/bx
441-8004002	4, Medium Adult	20/bx
441-8005002	5, Large Adult	20/bx



441-8005002




ICU Medical SACETT™ Endotracheal Tubes

The Portex® Suction Above the Cuff Endotracheal Tube (SACETT®) has been specifically designed to reduce VAP. The SACETT® combines the quality features of the Blue Line® range with the ability to remove secretions from above the cuff while leaving the tube in position.

- Posterior lumen with a subglottic evacuation opening above the cuff enables continuous aspiration of subglottic secretions
- The low pressure, high volume Soft-Seal® cuff is intended to minimize trauma to the trachea and protect the patient’s airway when inflated
- A spring-loaded, one-way valve helps ensure cuff inflation and integrity
- The large ergonomic pilot balloon makes it easier to distinguish between the sound of suction and that of a cuff leak

Item #	Description	ID	OD	Cuff Length	Qty
70-100/189/065	Murphy	6.5 mm	9.7 mm	23.0 mm	10/bx
70-100/189/070	Murphy	7.0 mm	10.4 mm	30.0 mm	10/bx
70-100/189/075	Murphy	7.5 mm	11.1 mm	30.0 mm	10/bx
70-100/189/080	Murphy	8.0 mm	11.9 mm	30.0 mm	10/bx
70-100/189/085	Murphy	8.5 mm	12.4 mm	30.0 mm	10/bx
70-100/189/090	Murphy	9.0 mm	12.8 mm	30.0 mm	10/bx
70-100/189/060	Murphy, Sub Glottic Suctioning	6.0 mm	—	—	10/bx




70-100/189/080



ICU Medical PORTEX® Clear PVC, Oral/Nasal, Soft Seal® Cuff Endotracheal Tubes

The PORTEX Soft Seal® Cuff combines the benefits of the profile cuff design with a larger cuff resting diameter to help offer the patient increased comfort over other products. The Soft Seal® cuff is made from soft PVC material, with the intention of reducing the risk of trauma. Strong, yet thin, cuffed material helps reduce the risk of folding and leakage.

Item #	ID	OD	Cuff Length	Qty
70-100/199/055EA	5.5 mm	7.6 mm	17 mm	20/bx
70-100/199/060EA	6.0 m	8.2 mm	23 mm	20/bx
70-100/199/065EA	6.5 mm	8.9 mm	23 mm	20/bx
70-100/199/070EA	7.0 mm	9.6 mm	30 mm	20/bx
70-100/199/075EA	7.5 mm	10.3 mm	30 mm	20/bx
70-100/199/080EA	8.0 mm	10.9 mm	30 mm	20/bx
70-100/199/085EA	8.5 mm	11.6 mm	30 mm	20/bx
70-100/199/090EA	9.0 mm	12.3 mm	30 mm	20/bx



70-100/199/080EA



Contact your local Tri-anim Account Manager for full range of ICU Medical Endotracheal Tubes.



our brand

Curaplex® Select Endotracheal Tubes with Stylet

- Standard 15 mm proximal connection with Murphy eye tip
- Sterile packaging for single-patient use
- Cuffs sized proportionally to the tube size
- Latex-free

Uncuffed Tube			Cuffed Tube		
Item #	ID	Qty	Item #	ID	Qty
301-2113-10325	2.5 mm	10/bx	301-2113-10225	2.5 mm	10/bx
301-2113-10330	3.0 mm	10/bx	301-2113-10230	3.0 mm	10/bx
301-2113-10335	3.5 mm	10/bx	301-2113-10235	3.5 mm	10/bx
301-2113-10340	4.0 mm	10/bx	301-2113-10240	4.0 mm	10/bx
301-2113-10345	4.5 mm	10/bx	301-2113-10245	4.5 mm	10/bx
301-2113-10350	5.0 mm	10/bx	301-2113-10250	5.0 mm	10/bx
			301-2113-10255	5.5 mm	10/bx
			301-2113-10260	6.0 mm	10/bx
			301-2113-10265	6.5 mm	10/bx
			301-2113-10270	7.0 mm	10/bx
			301-2113-10275	7.5 mm	10/bx
			301-2113-10280	8.0 mm	10/bx
			301-2113-10285	8.5 mm	10/bx
			301-2113-10290	9.0 mm	10/bx
			301-2113-10200	10.0 mm	10/bx



301-2113-10280



AirLife® Endotracheal Tubes

- Oral and nasal use
- Clear, flexible PVC
- Smooth Murphy Eye for alternative ventilation port
- Standard 15mm ISO connector included with each tube
- Disposable, latex-free, single-use, sterile

Uncuffed Tube					Cuffed Tube				
Item #	Fr	ID	Length	Qty	Item #	Fr	ID	Length	Qty
792-1-7330-20EA	8	2.0 mm	145 mm	10/bx	792-1-7333-30	12	3.0 mm	165 mm	10/bx
792-1-7330-25EA	10	2.5 mm	150 mm	10/bx	792-1-7333-35	14	3.5 mm	188 mm	10/bx
792-1-7330-30EA	12	3.0 mm	165 mm	10/bx	792-1-7333-40EA	16	4.0 mm	210 mm	10/bx
792-1-7330-35EA	14	3.5 mm	185 mm	10/bx	792-1-7333-45EA	18	4.5 mm	225 mm	10/bx
792-1-7330-40EA	16	4.0 mm	210 mm	10/bx	792-1-7333-50EA	20	5.0 mm	245 mm	10/bx
792-1-7330-45EA	18	4.5 mm	225 mm	10/bx	792-1-7333-55EA	22	5.5 mm	275 mm	10/bx
792-1-7330-50EA	20	5.0 mm	245 mm	10/bx	792-1-7333-60EA	24	6.0 mm	285 mm	10/bx
792-1-7330-55EA	22	5.5 mm	275 mm	10/bx	792-1-7333-65EA	26	6.5 mm	295 mm	10/bx
					792-1-7333-70EA	28	7.0 mm	310 mm	10/bx
					792-1-7333-75EA	30	7.5 mm	315 mm	10/bx
					792-1-7333-80EA	32	8.0 mm	330 mm	10/bx
					792-1-7333-85EA	34	8.5 mm	330 mm	10/bx
					792-1-7333-90EA	36	9.0 mm	330 mm	10/bx
					792-1-7333-95EA	38	9.5 mm	330 mm	10/bx
					792-1-7333-10EA	40	10.0 mm	330 mm	10/bx



792-1-7330-55



AirLife® Cuffed Endotracheal Tube with Stylet

- Oral and nasal use
- Clear, flexible PVC
- Smooth Murphy Eye for alternative ventilation port
- High volume, low pressure cuff
- Large barrel cuff
- Standard 15mm ISO connector included with each tube
- Disposable, latex-free, single-use, sterile

Item #	Fr	ID	Length	Qty
792-1-7343-30	12	3.0 mm	165 mm	10/bx
792-1-7343-35	14	3.5 mm	185 mm	10/bx
792-1-7343-40	16	4.0 mm	210 mm	10/bx
792-1-7343-45	18	4.5 mm	225 mm	10/bx
792-1-7343-50EA	20	5.0 mm	245 mm	10/bx
792-1-7343-55EA	22	5.5 mm	275 mm	10/bx
792-1-7343-60EA	24	6.0 mm	285 mm	10/bx
792-1-7343-65EA	26	6.5 mm	295 mm	10/bx
792-1-7343-70EA	28	7.0 mm	310 mm	10/bx
792-1-7343-75EA	30	7.5 mm	315 mm	10/bx
792-1-7343-80EA	32	8.0 mm	330 mm	10/bx
792-1-7343-85EA	34	8.5 mm	330 mm	10/bx
792-1-7343-90EA	36	9.0 mm	330 mm	10/bx



792-1-7343-75EA



AirLife® Nasal Pre-formed Endotracheal Tubes

- Designed to direct tube over patient's forehead, reducing pressure on nares
- Allows circuit to be positioned out of surgical field; ideal for oral and maxillofacial surgery applications
- 15mm male fitting included
- Smooth beveled tip for atraumatic introduction
- Disposable, latex-free, single-use, sterile



792-1-7323-55

Uncuffed Tube					Cuffed Tube				
Item #	Fr	ID	Length	Qty	Item #	Fr	ID	Length	Qty
792-1-7322-30	12	3.0 mm	160 mm	10/bx	792-1-7323-40	16	4.0 mm	200 mm	10/bx
792-1-7322-35	14	3.5 mm	180 mm	10/bx	792-1-7323-45	18	4.5 mm	220 mm	10/bx
792-1-7322-40	16	4.0 mm	200 mm	10/bx	792-1-7323-50	20	5.0 mm	240 mm	10/bx
792-1-7322-45	18	4.5 mm	220 mm	10/bx	792-1-7323-55	22	5.5 mm	270 mm	10/bx
792-1-7322-50	20	5.0 mm	240 mm	10/bx	792-1-7323-60	24	6.0 mm	280 mm	10/bx
792-1-7322-55	22	5.5 mm	270 mm	10/bx	792-1-7323-65	26	6.5 mm	290 mm	10/bx
					792-1-7323-70	28	7.0 mm	300 mm	10/bx
					792-1-7323-75	30	7.5 mm	310 mm	10/bx
					792-1-7323-80	32	8.0 mm	320 mm	10/bx



AirLife offers comprehensive portfolio of Endotracheal Tubes. For more information regarding additional sizes and styles, contact your local Tri-anim Account Manager.




ACCUCUFF™ CUFF PRESSURE INDICATOR

Representing a significant advance in cuff pressure monitoring, the AccuCuff™ cuff pressure indicator is one of our most innovative airway safety devices.



Calibrated to 20-29 cm H2O



Lightweight, efficient design



our brand


Curaplex® AccuCuff™ Cuff Pressure Indicator

- Continuous visual display eliminates the need for spot checks
- Negative pressure graduation ensures cuff is fully deflated, reducing the risk of trauma during removal
- Lightweight device reduces drag and potential tube displacement
- Suitable for use on a wide range of airway devices extending usage into many clinical areas
- Simple design for ease of use
- Pressure can be quickly and easily adjusted avoiding excessive cuff inflation
- Individually wrapped
- MRI compatible
- Single use
- Latex-free
- 5 year shelf life allows long term storage and reduces waste

Item #	Description	Qty
301-AC0100B-BX	Curaplex® AccuCuff™ ETT/Tracheal Cuff Pressure Indicator, Luer Lock (Blue)	50/bx
301-670249-KIT	Curaplex® AccuCuff™ ETT/Tracheal Cuff Pressure Indicator, Luer Lock (Blue)	5ea/kit
301-670252-KIT	Curaplex® AccuCuff™ ETT/Tracheal Cuff Pressure Indicator, Luer Lock (Blue) & Slip Tip Syringe KIT	10ea/kit



ETT/Tracheal Cuff Pressure Indicator 20-29cm H₂O



Intersurgical® InterForm Stylet

InterForm is a malleable endotracheal stylet, allowing the user to form the endotracheal tube into a suitable shape to ease insertion and give more control of the airway device.

- Aluminium core allows user to easily form suitable shapes
- Low-friction surface for easy insertion and removal
- Single-use, latex-free and supplied sterile in easy-to-open individual packs

Item #	Fr	Diameter	Size	Qty
441-8080006	6	2.0 mm	275 mm	10/bx
441-8080010	10	3.3 mm	340 mm	10/bx
441-8080014	14	4.7 mm	340 mm	10/bx

AirLife® Malleable Aluminum Stylet

- Smooth, rounded end
- Satin finish to reduce friction
- Individually wrapped
- Disposable, latex-free and sterile

Item #	Description	Fr	Size	Qty
792-9-0204-11	Pediatric Stylet	6	11"	25/bx
792-9-0204-25	Adult Stylet	9	16"	25/bx

AirLife® Aluminum Stylet

- Aluminum
- Disposable, latex-free

Item #	Description	Fr	Size	Qty
792-9-0204-00	Adult Stylet	10	16"	1/ea

AirLife® SunSlip Coated Stylet

- Sealed in a PVC sheath
- Malleable aluminum
- Disposable, latex-free and sterile

Item #	Description	Fr	Size	Qty
792-9-0203-06	Infant Stylet	6	10"	20/bx
792-9-0203-10	Child Stylet	10	11-1/2"	20/bx
792-9-0203-14	Adult Stylet	14	12"	20/bx

Intersurgical® InterGuide Bougie Introducer

InterGuide is a flexible tracheal tube introducer that allows positive location of the trachea and subsequent placement of the endotracheal tube in difficult airway situations. InterGuide is single-use, latex-free and supplied sterile in easy to open individual packs.

Item #	Tip	Fr	Diameter	Size	Qty
441-8070006	Coudé Tip	6	2.0 mm	530 mm	10/bx
441-8070010	Coudé Tip	10	3.3 mm	700 mm	10/bx
441-8070015	Coudé Tip	15	5.0 mm	700 mm	10/bx

AirLife® Bougie-To-Go™ ET Introducer

- Unwraps without distorting introducer shape; keeps tip in correct insertion position
- Low density polyethylene provides proper stiffness for ease of insertion
- Disposable, latex-free and sterile

Item #	Description	Tip	Fr	Size	Qty
792-9-0212-82	Adult	Coudé Tip	15	60 cm	10/bx

AirLife® Original Blue Bougie Introducer


- Provides proper stiffness for ease of insertion
- Calibrated for easy observation of insertion
- Coudé tip facilitates insertion in adults
- Pediatric size has calibration marks on both ends allowing use as straight or Coudé tip introducer
- Individually-wrapped
- Disposable, latex-free and sterile

Item #	Description	Tip	Fr	Size	Qty
792-9-0211-70	Pediatric	Coudé/ Straight Tip	10	70 cm	10/bx
792-9-0212-72	Adult	Straight Tip	15	70 cm	10/bx
792-9-0212-70	Adult	Coudé Tip	15	70 cm	10/bx

AirLife® Malleable Bougie Introducer

- Combines benefits of an introducer and stylet into one unit
- Internal wire provides user the ability to pre-bend the stylet
- Graduation marks every 5cm
- Can be used with ET tubes, rigid bronchoscope, laryngeal type masks and double lumen endobronchial tubes
- Disposable, latex-free and sterile


Item #	Description	Fr	Size	Qty
792-9-0213-82	Pediatric	10	70 cm	10/bx
792-9-0213-92	Adult	15	70 cm	10/bx



792-9-0212-82



792-9-0211-70



792-9-0212-72

792-9-0212-70

792-9-0203-14



792-9-0203-06



AirLife® Adjustable Flange Nasopharyngeal Airways

- Flexible neoprene performs like latex airways
- Soft material prevents damage to nasal passage
- Soft, adjustable flange
- Blunt tip
- Individually wrapped
- Disposable, latex-free, single-use, sterile


Item #	Fr	D	Qty
792-1-5072-12EA	12	3 mm	10/bx
792-1-5072-14EA	14	3.5 mm	10/bx
792-1-5072-16EA	16	4 mm	10/bx
792-1-5072-18EA	18	4.5 mm	10/bx
792-1-5072-20EA	20	5 mm	10/bx
792-1-5072-22EA	22	5.5 mm	10/bx
792-1-5072-24EA	24	6 mm	10/bx
792-1-5072-26EA	26	6.5 mm	10/bx
792-1-5072-28EA	28	7 mm	10/bx
792-1-5072-30EA	30	7.5 mm	10/bx
792-1-5072-32EA	32	8 mm	10/bx
792-1-5072-34EA	34	8.5 mm	10/bx
792-1-5072-36EA	36	9 mm	10/bx



AirLife® PVC Nasopharyngeal Airways

- Smooth surface for easy insertion
- Beveled tip
- 15mm fitting
- Individually wrapped
- Disposable, latex-free, single-use, sterile

Item #	Fr	D	Qty
792-1-5073-12	12	3 mm	10/bx
792-1-5073-14	14	3.5 mm	10/bx
792-1-5073-16	16	4 mm	10/bx
792-1-5073-18	18	4.5 mm	10/bx
792-1-5073-26	26	6.5 mm	10/bx
792-1-5073-28	28	7 mm	10/bx
792-1-5073-30	30	7.5 mm	10/bx
792-1-5073-32	32	8 mm	10/bx
792-1-5073-34	34	8.5 mm	10/bx



792-1-5073-32



AirLife® Latex-Free Trumpet Style (Robertazzi) Nasopharyngeal Airways

- Flexible neoprene performs like latex airways
- Soft material prevents damage to nasal passage
- Blunt tip
- Resists collapse or kinking
- Sets are ideal for crash carts
- Individually wrapped
- Disposable, latex-free, single-use, sterile

Item #	Description	Fr	D	Qty
792-1-5075-20EA	Airway, Robertazzi	20	5 mm	10/bx
792-1-5075-22EA	Airway, Robertazzi	22	5.5 mm	10/bx
792-1-5075-24EA	Airway, Robertazzi	24	6 mm	10/bx
792-1-5075-26EA	Airway, Robertazzi	26	6.5 mm	10/bx
792-1-5075-28EA	Airway, Robertazzi	28	7 mm	10/bx
792-1-5075-30EA	Airway, Robertazzi	30	7.5 mm	10/bx
792-1-5075-32EA	Airway, Robertazzi	32	8 mm	10/bx
792-1-5075-34EA	Airway, Robertazzi	34	8.5 mm	10/bx
792-1-5075-36EA	Airway, Robertazzi	36	9 mm	10/bx

792-1-5075-04	4 pk-sizes 20, 24, 28 & 32 plus Sterile Lub. Jelly
792-1-5075-06	6 pk-sizes 22, 24, 26, 28, 30 & 32 plus Sterile Lub. Jelly
792-1-5075-09	9 pk-sizes 20, 22, 24, 26, 28, 30, 32, 34 & 36 plus Sterile Lub. Jelly



our brand

Curaplex® Nasopharyngeal Airways

Curaplex® nasopharyngeal airways are latex-free, sterile, single-use and packaged individually. Each package features a color coded label to clearly denote French sizing. This unique new package design enables immediate and accurate airway selection.

Item #	Description	Fr	D	Qty
2021-14630	PVC	12	3.0 mm	10/bx
2021-14635	PVC	14	3.5 mm	10/bx
2021-14640	PVC	16	4.0 mm	10/bx
2021-14645	PVC	18	4.5 mm	10/bx
2021-14650	PVC	20	5.0 mm	10/bx
2021-14655	PVC	22	5.5 mm	10/bx
2021-14660	PVC	24	6.0 mm	10/bx
2021-14665	PVC	26	6.5 mm	10/bx
2021-14670	PVC	28	7.0 mm	10/bx
2021-14675	PVC	30	7.5 mm	10/bx
2021-14680	PVC	32	8.0 mm	10/bx
2021-14685-R	PVC	34	8.5 mm	10/bx
2021-14690	PVC	36	9.0 mm	10/bx

2021-18120	Neoprene	20	5.0 mm	10/bx
2021-18122	Neoprene	22	5.5 mm	10/bx
2021-18124	Neoprene	24	6.0 mm	10/bx
2021-18126	Neoprene	26	6.5 mm	10/bx
2021-18128	Neoprene	28	7.0 mm	10/bx
2021-18130	Neoprene	30	7.5 mm	10/bx
2021-18132	Neoprene	32	8.0 mm	10/bx
2021-18134	Neoprene	34	8.5 mm	10/bx
2021-18136	Neoprene	36	9.0 mm	10/bx



2021-14690

2021-18136



AirLife® PVC Guedel Airways

- Soft, smooth flexible clear PVC finish for patient safety
- Color-coded for easy size identification
- Rigid built-in bite block
- Individually wrapped
- Disposable, latex-free
- Non-sterile

Item # Pk / 10	Item # Box / 50	mm	Size	Color
792-1-5010-30	792-1-1521-30	30	000	Clear
792-1-5010-20	792-1-1521-20	40	00	Pink
792-1-5010-10	792-1-1521-10	50	0	Blue
792-1-5010-01	792-1-1521-01	60	1	Black
792-1-5010-02	792-1-1521-02	70	2	White
792-1-5010-03	792-1-1521-03	80	3	Green
792-1-5010-04	792-1-1521-04	90	4	Yellow
792-1-5010-05	792-1-1521-05	100	5	Red
792-1-5010-06	792-1-1521-06	110	6	Orange
792-1-5010-98	Complete PVC Guedel Airway Set, One of Each size			



our brand

Curaplex® Guedel Airways, Disposable

- The center passage permits suction catheters to clear secretions from the back of the throat
- Smooth, plastic finish allows simple insertion, minimizing tissue trauma
- Individually wrapped and color-coded for easy size identification
- Disposable for single patient use
- Non-sterile
- Serves multiple purposes as a biteblock and tongue depressor which avoids crimping and protects teeth

Item #	mm	Color	Qty
301-2010-35040	40	Pink	50/bg
301-2010-35050	50	Blue	50/bg
301-2010-35060	60	Black	50/bg
301-2010-35070	70	White	50/bg
301-2010-35080	80	Green	50/bg
301-2010-35090	90	Yellow	50/bg
301-2010-35100	100	Red	50/bg
301-2010-35110	110	Light Blue	50/bg



AirLife® Guedel Airway

- Firm, polyethylene plastic
- Rigid bite block
- Individually wrapped
- Color-coded
- Disposable, latex-free
- Non-sterile

Item # Pk / 10	Item # Box / 50	mm	Size	Color
792-1-1504-40EA	792-1-1505-40	40	00	Pink
792-1-1504-50EA	792-1-1505-50EA	50	0	Blue
792-1-1504-60EA	792-1-1505-60	60	1	Black
792-1-1504-70EA	792-1-1505-70	70	2	White
792-1-1504-80	792-1-1505-80EA	80	3	Green
792-1-1504-90	792-1-1505-90EA	90	4	Yellow
792-1-1504-99EA	792-1-1505-99	100	5	Red
792-1-1504-11	792-1-1505-11	110	6	Orange
792-1-1504-00	Complete Guedel Airway with Bite Block Set, One of Each size			




Intersurgical® Guedel Airways

One-piece Guedel airways combine two materials with very different physical characteristics in one moulding, resulting in the elimination of the separate bite block found in all conventional designs.

Item #	mm	Color	Qty
441-1000SD	35	Pink	50/ca
441-1100SD	50	Blue	50/ca
441-1110SD	55	Gray	50/ca
441-1111SD	65	White	50/ca
441-1111031	70	Yellow	50/ca
441-1112SD	80	Green	50/ca
441-1113SD	90	Orange	50/ca
441-1114SD	100	Red	50/ca
441-1115SD	120	Purple	50/ca



301-2010-34110



301-2010-34090

301-2010-34100

301-2010-34080



our brand

Curaplex® Berman Oral Airways, Disposable

- Holes through the center of the curve allow for additional air passage
- Open channels on each side to enable suctioning
- Flat surface for patient comfort
- Polypropylene material
- Color-coded for easy size identification
- Non-sterile, disposable for single patient use

Item #	mm	Color	Qty
301-2010-34040	40	Pink	50/bg
301-2010-34050	50	Blue	50/bg
301-2010-34060	60	Black	50/bg
301-2010-34070	70	White	50/bg
301-2010-34080	80	Green	50/bg
301-2010-34090	90	Yellow	50/bg
301-2010-34100	100	Red	50/bg
301-2010-34110	110	Light Blue	50/bg

792-1-1508-90EA



792-1-1508-11EA

792-1-1508-99EA



AirLife® Berman Airways

- Vented septum allows air flow from either side
- Firm, polyethylene plastic
- Color-coded for easy size identification
- Rigid built-in bite block
- Individually wrapped
- Disposable, latex-free, and non-sterile

Item # Pk / 10	Item # Box / 50	mm	Size	Color
792-1-1508-40EA	792-1-1506-40EA	40	00	Pink
792-1-1508-50EA	792-1-1506-50EA	50	0	Blue
792-1-1508-60EA	792-1-1506-60EA	60	1	Black
792-1-1508-70EA	792-1-1506-70EA	70	2	White
792-1-1508-80EA	792-1-1506-80EA	80	3	Green
792-1-1508-90EA	792-1-1506-90EA	90	4	Yellow
792-1-1508-99EA	792-1-1506-99EA	100	5	Red
792-1-1508-11EA	792-1-1506-11EA	110	6	Orange
792-1-1508-00	Complete Berman Airway Set, ten of each size			



AirLife® Williams Airway Intubator

- Use as an oropharyngeal airway to intubate the trachea
- Becomes a guide for fiber optic bronchoscope placement
- Rigid built-in bite block
- Disposable, latex-free, and non-sterile

Item #	Description	Size	Qty
792-1-1503-09	Airway Intubator for Adult Females	9 cm	10/bx
792-1-1503-10	Airway Intubator for Adult Males	10 cm	10/bx

792-1-1503-09



792-1-1503-10

AirLife® SunStim™ Peripheral Nerve Stimulators

SunStim™:

- Stand-by position: no stimulus pulses are generated
- Twitch: 1 pulse per second
- Tetanus: 100 Hz (100 pulses/s)
- Built-in train of four: 4 pulses in a period of 2 seconds

SunStim™ Plus:

- Double burst: 2 60-ms bursts of 50 Hz separated by .75 s
- Twitch: 1 pulse per second
- Tetanus: 50 Hz (50 pulses/s) or 100 Hz (100 pulses/s)
- Built-in train of four: 4 pulses in a period of 2 seconds

Item #	Description	Qty
792-8-1053-60	SunStim Peripheral Nerve Stimulator	1/ea
792-8-1053-62	SunStim Plus Peripheral Nerve Stimulator	1/ea
792-8-1053-76	Ball Electrode and Hollow Pin	1/ea
792-8-1051-78	6' Lead, Snap and Shielded Pin	1/ea
792-8-1051-79	6' Lead, Clip and Shielded Pin	1/ea



792-8-1053-60



792-8-1053-62

AirLife® SunStim™ Pro Peripheral Nerve Stimulator

Ergonomically designed stimulator that delivers reliable, adjustable current with a digital screen for quick visual assessment of functionality. Features a removable silicone wrap which enhances grip and helps to protect the device. The silicone wrap comes in a variety of colors.

SunStim™ Pro Functionality:

- Double burst: 2 60-ms bursts of 50 Hz separated by 0.75 s
- Twitch: 1 pulse per second
- Tetanus: 50 Hz (50 pulses/s) or 100 Hz (100 pulses/s)
- Built-in train of four: 4 pulses in a period of 2 seconds
- Stimulus pulse: square wave monophasic pulses (200-μs duration)

Item #	Description	Qty
792-8-1053-63	SunStim™ Pro Peripheral Nerve Stimulator	1/ea
792-8-1053-76	Ball Electrode and Hollow Pin	1/ea
792-8-1051-78	6' Lead, Snap and Shielded Pin	1/ea
792-8-1051-79	6' Lead, Clip and Shielded Pin	1/ea
792-8-1053-95	Pole Mount Bracket for 8-1053-63	1/ea
792-8-1053-87	Zippered Pouch 8.5" x 6"	1/ea
792-8-1053-GR	Silicone Pro Wrap Gray	1/ea
792-8-1053-BL	Silicone Pro Wrap Blue	1/ea
792-8-1053-PL	Silicone Pro Wrap Purple	1/ea
792-8-1053-FS	Silicone Pro Wrap Fuchsia	1/ea
792-8-1053-PK	Silicone Pro Wrap Pink	1/ea
792-8-1053-RD	Silicone Pro Wrap Red	1/ea
792-8-1053-OR	Silicone Pro Wrap Orange	1/ea
792-8-1053-GN	Silicone Pro Wrap Green	1/ea
792-8-1053-TL	Silicone Pro Wrap Teal	1/ea




792-8-1053-63



792-8-1053-GN



792-8-1053-87



792-8-1053-95



715-PM5300

Custom Hose Configurations

Precision Medical offers a wide variety of custom hose configurations to meet all OR needs, such as WAGD, along with the appropriate fittings.

<http://hoseconfig.precisionmedical.com>



Precision Medical High Performance Blenders

Precision Medical offers high performance blenders that maintain FiO₂ even at low flows, thereby reducing the risk of OR fires. These user-friendly units are durable, versatile and offer cost saving opportunities through reduced O₂ use during procedures. Made in USA.

Item #	Description	Qty
715-PM5200	High Flow Blender	1/ea
715-PM5200M	High Flow Blender with Manifold	1/ea
715-PM5201	High Flow Blender w/ outlets on L/R sides	1/ea
715-PM5300	Low Flow Blender	1/ea
715-PM5300M	Low Flow Blender with Manifold	1/ea



For a comprehensive offering of Precision Medical's blenders, suction, flowmeters, medical fittings, custom hose configurations and more, contact your local Tri-anim Account Manager.



Reduce the Risk of OR Fires

Procedures involving the head, neck and upper chest/above T5 are at greatest risk due to use of an ignition source in proximity to an oxidizer. The APSF recommends delivery devices, such as blenders or common gas outlets, maintain oxygen below 30% to safely manage ignition sources.



Maxtec® MaxBlend™ 2 Air-Oxygen Blender with Built-in Monitor

MaxBlend™ 2, which was designed with years of customer feedback in mind. This easy-to-use and low-maintenance blender has a host of new features and improvements, including a backlit LCD and acrylic flow meter, smart alarms, a DC power port, extended battery life, and more which are all backed by the Maxtec three-year warranty.

- Three-in-one blender
- Helps eliminate cross-contamination
- Provides measurements free from temperature, pressure, and humidity
- 15% to 99% low alarm capabilities
- Multiple flow meter configurations and bleed control
- Smart alarms – sets alarms 3% +/- current reading

Item #	Description	LPM	Qty
203-R229P02-001	MaxBlend 2 Blender	0-70	1/ea
203-R229P01-010	MaxBlend 2 Blender	0-15	1/ea
203-R229P01-011	MaxBlend 2 Blender	0-3	1/ea



For a comprehensive offering of Maxtec's blenders, flowmeters, analyzers and accessories, please contact your local Tri-anim Account Manager.



INTERSURGICAL® CO₂ ABSORBENTS

Intersurgical offers a complete range of medical grade CO₂ absorbents and standard soda lima, that meet all clinical requirements and safety needs within anesthesia.



What is in the products?

	Spherasorb™	LoFloSorb™	Intersorb Plus™	Other conventional soda limes
Calcium Hydroxide	93.5%	92.5%	97%	97%
Sodium Hydroxide	1.5%	Nil	3%	3%
Silica Synthetic/ amorphous	Nil	7.5%	Nil	Nil
Zeolite	5%	Nil	Nil	Nil
Indicator	0.03%	0.03%	0.03%	0.03%



Intersurgical® CO₂ Absorbents

Item #	Description	Size	Qty
441-2186030	Spherasorb™ Drum, Pre-filled Cartridge	1 kg	1/ca
441-2191	Spherasorb™ Pyramid, Pre-filled, Disposable	1 kg	6/ca
441-2188030	LoFloSorb™ Drum, Pre-filled Cartridge	1 kg	10/ca
441-2193	LoFloSorb™ Pyramid, Pre-filled, Disposable	1 kg	6/ca
441-2196	Spherasorb™ IS Can disposable CO ₂ absorber	1 kg	6/ca
441-2198	LoFloSorb™ IS Can disposable CO ₂ absorber	1 kg	6/ca
441-2175030	Spherasorb™ Jerican	5 L	2/ca

Tracheostomy T-Pieces

7.6mm Port
22 OD
22 OD/15 ID
22 OD

441-1984000

22 OD
22 OD
22 OD/15 ID

441-1980000

Monitoring Connectors

Luer Port & Cap
15 OD
22 OD/15 ID

441-2714000

Female Luer Fitting
7.6mm Stem
441-27106000

Luer Port & Cap
15 OD
22 OD/15 ID

441-2713000

MDI Connectors

MDI Insert
22 OD
22 OD

441-1966000

MDI Insert
22 OD
22 ID

441-1964001

MDI Insert
15 OD
15 ID

441-1947000

MDI Insert
15 OD
15 ID

441-1947003

Elastomeric Connectors

For 6-9mm Spigot
15 ID
441-1702000

For 9-11mm Spigot
15 ID
441-17036000

For 22mm lipped or non ISO connectors
22 ID
441-1704000

For 22mm lipped or non ISO connectors
22 ID
441-1706

For 22mm lipped or non ISO connectors
22 ID
441-1705

For 22mm lipped or non ISO connectors
22 ID
441-1701000

Scavenging Connectors

22 OD
19 OD
441-1975000

22 OD
30 OD
441-1970000

22 OD
30 ID
441-1971000

Swivel Elbows

Flip top cap with 7.6mm port
15 OD Swivel
22 OD/15 ID Swivel
441-1899000

15 OD Swivel
22 OD/15 ID Swivel
441-1897000

Double Flip top cap with seal
15 OD Swivel
22 OD/15 ID Swivel
441-1898000-30

Oxygen Connectors & Adaptors

22 OD
6mm Oxygen Stem
441-1968000

6mm Oxygen Stem
22 ID
441-1568000

6.5mm Oxygen Stem
15 ID/22 OD
441-1974000

6mm swivel system
22 OD
22 ID
441-1963000

Venturi Valves

22 ID
6mm Stem
441-0024030

22 ID
6mm Stem
441-0028030

22 ID
6mm Stem
441-0031030

22 ID
6mm Stem
441-0035030

22 ID
6mm Stem
441-0040030

22 ID
6mm Stem
441-0060030

22 ID
6mm Stem
441-0070030

Y-Pieces

22 OD
22 OD
22 OD/15 ID
441-19000000

22 OD
7.6mm port
6mm stem
22 OD/15 ID
441-1901000

22 OD Swivel
Luer Lock Port
22 OD/15 ID Swivel
22 OD Swivel
441-19296000

15 OD
15 OD
441-1902000

15 OD
7.6mm ports
22 OD/15 ID
441-1903000

22 OD
7.6mm ports
22 OD/15 ID
441-1927000

22 OD
22 OD
22 OD/15 ID
441-19266000

15mm Straight Connectors

15 OD
15 OD
441-1945000

15 OD
15 ID
441-1946000

7.6mm Port
15 OD
15 ID
441-1947000

6mm Swivel Stem
15 OD
15 ID
441-1948030

15 ID
15 ID
441-1949000

15 OD
22 OD
441-1943000

15 OD
22 ID
441-1962000

15 ID
22 ID
441-1979000

Luer Port and Cap
15 OD
15 ID
441-2712030

Fixed Elbows

15 OD
22 OD/15 ID
441-1997000-30

Flip top cap with 7.6mm port
15 OD
22 OD/15 ID
441-1993000

7.6mm Port
22 OD
22 OD/15 ID
441-1995000

6mm Stem
22 OD
22 OD/15 ID
441-1990000

7.6mm Port
22 ID
22 OD/15 ID
441-1893

22 ID
22 OD
441-1992000

6mm Stem
22 ID
22 OD
441-1991000

7.6mm Port
22 ID
22 ID
441-1891000

22mm Straight Connectors

22 OD
22 OD
441-1969000

22 OD
22 ID
441-1961000

7.6mm Port
22 OD
22 OD
441-1966000

6mm Swivel Stem
22 OD
22 OD
441-1965000

7.6mm port
22 OD
22 ID
441-1964000

22 OD
22 OD/15 ID
441-1969000

22 ID
22 ID
441-1967000

22 ID
441-1978000

Mouthpiece
22 OD
22 OD
441-1510000

Mouthpiece
22 OD
441-1930030

Mouthpiece
22 ID
441-1931030

22 OD/15 ID
22 OD/15 ID
441-1933000

Tube Management Solutions

Tube management clips
441-2316000

Notch feature
22 mm
441-2316000

Notch feature
22 mm
441-1978000

T-Pieces

15 OD
15 ID
15 ID
441-1802

15 OD
15 ID
22 OD
441-18006000

15 OD
22 OD
15 OD
441-1988000

15 OD
22 OD
22 OD
441-1989000

22 OD
15 ID
22 OD
441-1959000

22 OD/15 ID
22 ID
22 ID
441-1983000

22 OD
22 ID
22 ID
441-1981000

22 OD
22 OD/15 ID
22 ID
441-1985000

22 OD/15 ID
22 OD
22 ID
441-1986

22 OD
22 ID
22 ID
441-1982000

Superset Tube Catheter Mounts

22 ID
22 OD/15 ID
441-3504000 | STERILE OPTION: 441-3504000S

22 ID
22 OD/15 ID
441-3502000

22 ID
Swivel
Flip top cap with 7.6mm port
22 OD/15 ID
441-3520000

22 ID
Swivel
Double Flip top cap with seal. Covers 7.6mm and 9.5mm applications
22 OD/15 ID
441-3521000

22 ID
Swivel
22 OD/15 ID
441-3512EA

Luer port and cap
22 ID
22 OD/15 ID
441-3514031

22 ID
15 ID
441-9123

15 OD
7.6mm Port
15 ID
441-3532000

Flexible Tube Catheter Mounts

22 ID
22 OD/15 ID
441-3511

22 ID
Swivel
22 OD/15 ID Swivel
441-3506

Directional Valves

22 OD
22 OD
441-1950000

22 OD
22 ID
441-1951

22 OD
22 OD
22 OD
441-1953

22 OD
22 OD
22 OD/15 ID
441-1954

22 OD
22 ID
22 OD
441-1952000

22 OD
22 OD
22 OD/15 ID
441-1958



FILTERS AND HMEs

Curaplex® offers a comprehensive portfolio of filtration products and heat and moisture exchangers (HME).



Curaplex® Filters and HMEs

	FILTERS			HME				HMEF		
Item #	301-439	301-431	301-435	301-1014	301-422	301-427	301-429	301-437	301-424	301-423
Description	Aero-Pro™ HEPA Light with Port	Aero-Pro™ HEPA Light Machine	BFV/VFE Depth Filter	ClearPath Aerosol HME	Aero-Sat™ 2 Straight	Aero-Sat™ 1 Straight	Aero-Sat™ Mini	Aero-Sat™ HEPA	Aero-Sat™ Filter Compact Straight w/ Flex Tube	Aero-Sat™ Filter Compact Straight
Moisture Output (mg/L @ TV(mL))	—	—	—	34.2 @ 500 31.3 @ 1000	31 @ 500	28 @ 250	18 @ 20	30 @ 500	30 @ 1000	30 @ 1000
Flow Resistance (cmH ₂ O @ L/min)	2.32 @ 60	2.0 @ 60	2.1 @ 60	0.63 @ 30	1.15 @ 60	.41 @ 20	1.05 @ 10	2.48 @ 60	2.45 @ 60	2.45 @ 60
Dead Space (mL)	103.8	80	29	80.4	31.0	13.2	3.1	125.6	46.2	46.2
Weight (g)	40.1	37	19.8	34	22.5	9.8	4.1	49.6	33.2 (37.9 W/TUBE)	33.2
Tidal Volume (mL)	300-1200	300-1200	>125	250-1000	150-1500	50-600	15-50	300-1200	150-1000	150-1000
Media	Hydrophobic Paper	Hydrophobic Paper	Electrostatic Filter	Foam	HME Paper	HME Paper	HME Paper	HME Paper/ Hydrophobic Paper	HME Paper/ Electrostatic Filter	HME Paper/ Electrostatic Filter
Port	Yes	No	No	—	No	No	No	Yes	Yes	Yes
Bacterial Efficiency (%)	99.99999	99.99999	99.999	—	—	—	—	99.99999	99.999	99.999
Viral Efficiency (%)	99.9999	99.9999	99.999	—	—	—	—	99.9999	99.99	99.99



ICU Medical Portex® HME Filter, HEPA Straight

ICU Medical Portex® HMEs provide efficient heat and moisture recovery. All HMEs exhibit low resistance to flow and high moisture output which helps to minimize the loss of heat and moisture in ventilated or spontaneously breathing patients. All HMEs and filters have ISO standard tapered ends to assure a consistent leak-free seal.

Item #	70-002873
Description	Heat Moisture Exchanger Filter HMEF, HEPA, Pleated Paper, Tidal Volume 200-1500ml, Straight
Tidal Volume /ml	200-1500
Resistance /cm H ₂ O @ L/min	1.5 @ 30
Humidity Output /mg H ₂ O/L @ stated Vt	25 @ 500
Deadspace /ml 76 52 55	52
Weight /g	40
Medium	Pleated Paper
Filter	Hydrophobic
Filter Efficiency	99.9999%
Units/Case	50/ca



70-002873



ICU Medical Portex® HME Filter, Mechanical HEPA with Gas Port, Straight

Item #	70-002874
Description	Heat Moisture Exchanger Filter HMEF, Mechanical HEPA, Pleated Paper, HME Foam, Gas Port, Straight
Tidal Volume /ml	150-1500
Resistance /cm H ₂ O @ L/min	1.4 @ 30
Humidity Output /mg H ₂ O/L @ stated Vt	31 @ 500
Deadspace /ml 76 52 55	55
Weight /g	40
Medium	Pleated Paper, HME Foam
Filter	Hydrophobic
Filter Efficiency	99.9999%
Units/Case	50/ca



70-002874



AirLife® Hygroscopic Condenser Humidifier (HCH)

Hygroscopic condenser humidifiers (also called heat and moisture exchangers - HME/HMEF) effectively provide heat and moisture delivery/recovery for ventilated and spontaneously breathing patients. Available in a variety of sizes and shapes, HCHs are compatible with leading breathing systems and come with or without filtration material. Latex-free. Single-patient use. Disposable.

With Filter (HMEF)		
Item #	Description	Qty
87-FH603005EA	Round, Large, without Port	50/ca
87-FH603008	Round, Large, with Port	50/ca
87-FH603009EA	Swivel, with Port	50/ca
Without Filter (HME)		
87-FH603021	Straight, Small, with Port	50/ca
87-FH603023EA	Straight, Small, without Port	50/ca

 For a complete listing of AirLife HME/HMEFs, contact your Tri-anim Account Manager.

Intersurgical® Filta-Therm™ Heat and Moisture Exchanging Filters (HMEF)

The range of Heat and Moisture Exchanging Filters (HMEFs) combine the filtration efficiency of dedicated breathing filters with optimum moisture return.

Item #	Description	Qty
441-1942030	Filta-Therm	40/ca
441-1941030	Filta-Therm with luer lock port	40/ca
441-1942000	Filta-Therm with flextube	20/ca
441-19416500	Filta-Therm with luer lock port, flextube	20/ca



441-1942030

AirLife® Breathing Bags

- Single-use, latex-free
- Pleated bag design

Item #	Size	Qty
87-HM025	¼ Liter Bag	20/ca
87-HM05	½ Liter Bag	20/ca
87-HM10EA	1 Liter Bag	20/ca
87-HM20EA	2 Liter Bag	20/ca
87-HM30EA	3 Liter Bag	20/ca



AirLife® Flat Breathing Bags

- Reusable, latex-free
- Fabricated from the finest neoprene
- 22 mm bushing

Item #	Size	Qty
792-4-1041-05	1/2 Liter Bag	1/ea
792-4-1041-10	1 Liter Bag	1/ea
792-4-1041-15	1-1/2 Liter Bag	1/ea
792-4-1041-20	2 Liter Bag	1/ea
792-4-1041-30	3 Liter Bag	1/ea
792-4-1041-40	4 Liter Bag	1/ea
792-4-1041-50	5 Liter Bag	1/ea



AirLife® Filters

Item #	Description	Qty
963-6218	HME Filter, Angled Inlet with CO ₂ Port	50/ca
963-6219	In-Line Mini HME with CO ₂ Port	50/ca
963-6220	In-Line Mini HMEF with CO ₂ Port	50/ca
963-6221	In-Line Maxi HME with CO ₂ Port	50/ca
963-6229	HME, Angled Inlet, with CO ₂ Port	50/ca
963-6232	Infant, HMEF with CO ₂ Port, 15 OD x 15 mm ID	50/ca
963-6379	In-Line Maxi HMEF, No Port	50/ca



Intersurgical® Filta-Guard™ and Clear-Guard™

The Intersurgical range of breathing filters have been designed for the protection of the patient, breathing system and equipment.

Item #	Description	Qty
441-1944030	Filta-Guard	40/ca
441-1544030	Clear-Guard 3 filter + luer lock port	40/ca
441-1644030	Clear-Guard Midi + luer lock port	50/ca

 For a complete listing of Intersurgical HME/HMEFs, contact your Tri-anim Account Manager.





441-2158009

Intersurgical® Multipac Anesthetic Breathing Circuits

Cost savings are a major factor in the modern surgical environment. This not only affects products, but also storage requirements and anesthesia set-up time. Intersurgical Incorporated offers a completely integrated anesthesia system to address all of these concerns. Multipac is a one-a-day breathing circuit that may be used by up to six patients. Each breathing system contains two patient packs - see contents below. By selecting one of our eight Multipac configurations, most circuit needs can be met based on your specifications.

Item #	Breathing Circuit Type			Patient Pack Contents (Contains 2 Packs)			Breathing System Configuration						Qty
	Compact II	UniFlow™	Flextube™	ClearLite Mask (Size 4)	Economy Mask (Size 4)	Filta-Therm HMEF & Elbow	Length	Parallel Y-Piece	End Cuffs	3 L Reservoir Bag (latex free)	Monitoring Line Length	Integral Monitoring Line	
441-2158009	✓			✓		✓	118 inches	✓	✓	✓	3 m		10/bx
441-2158010	✓				✓	✓	118 inches	✓	✓	✓	3 m		10/bx
441-2900077		✓		✓		✓	63 inches		✓	✓	3 m		10/bx
441-2900078		✓			✓	✓	63 inches		✓	✓	3 m		10/bx
441-2902022		✓		✓		✓	63 inches		✓	✓	.9 m	✓	10/bx
441-2902023		✓			✓	✓	63 inches		✓	✓	.9 m	✓	10/bx
441-2010235			✓	✓		✓	63 inches	✓	✓	✓	3 m		10/bx
441-2010261			✓		✓	✓	63 inches	✓	✓	✓	3 m		10/bx



441-72936000

Intersurgical® Multipac Patient Packs

Additional options are available by selecting one of our eight Patient Packs that will include our high efficiency Filta- Therm™ HMEF, an anesthesia mask elbow and your choice of anesthetic face masks: ClearLite™ or Economy anesthetic face mask. Packs must be exchanged between each patient.


By using Multipac systems, cost savings can be realized through the product, the reduction of storage needs and a reduction in waste costs.

Item #	Description	Size	Qty
441-72936000	Filta-Therm HMEF, anesthesia mask elbow and ClearLite anesthetic face mask	3	20/bx
441-72946000	Filta-Therm HMEF, anesthesia mask elbow and ClearLite anesthetic face mask	4	20/bx
441-72956000	Filta-Therm HMEF, anesthesia mask elbow and ClearLite anesthetic face mask	5	20/bx
441-72966000	Filta-Therm HMEF, anesthesia mask elbow and ClearLite anesthetic face mask	6	20/bx
441-15146000	Filta-Therm HMEF, anesthesia mask elbow and Economy anesthetic face mask	3	20/bx
441-15156000	Filta-Therm HMEF, anesthesia mask elbow and Economy anesthetic face mask	4	20/bx
441-15166000	Filta-Therm HMEF, anesthesia mask elbow and Economy anesthetic face mask	5	20/bx
441-15176000	Filta-Therm HMEF, anesthesia mask elbow and Economy anesthetic face mask	6	20/bx

Intersurgical® UniFlow™ Coaxial Breathing Systems

UniFlow™ is a convenient tube-within-a-tube breathing system designed with low flow anesthesia in mind, offering warmer, more humid inspiratory gases, and the option of an integrated respiratory gases monitoring line. UniFlow can now be used on pediatrics to as low as 5 kg in weight. Each UniFlow system is supplied complete with a leak tester.


Item #	ID/OD	L	Qty
441-29006000	30 mm	1.6 M with 0.5 M Limb	10/ca
441-29016000	30 mm	2.4 M with 0.5 M Limb	10/ca
441-29026000	30 mm with Integrated Monitor Line	1.6 M with 0.5 M Limb	10/ca
441-29036000	30 mm with Integral Monitoring Line	2.4 M	10/ca



AirLife® Anesthesia Breathing Circuits

High quality tubing with exceptional flexibility and drapability. ISO-standard connections for secure, leak resistant connections. Available in a wide variety of lengths and configurations.

Item #	Description	Size	Bag	Qty
963-8533	Dual-Limb, 60 in, Parallel Wye, Gas Sampling Elbow an Line, Filter	Adult	3 L	10/ea
963-7-SP5651	Dual-Limb, 72 in, Bifurcated Wye, Elbow w/ Red Cap, Water Trap, 22 mm Tubing	Adult	—	10/ea
963-7-8113-10-68	Dual-Limb, 96 in, Parallel Wye, Gas Sampling Elbow, Filter, 3 L Bag	Adult	3 L	15/ea
963-8090	Dual-Limb, 120in, Parallel Wye, Angled HME w/ Port	Pediatric	1 L	20/ea
963-7-89-11-10	Dual-Limb, 36 in, Gas Sampling Elbow, Parallel Wye	Pediatric	1 L	30/ea
963-9096	UniLim, 120 in, Gas Sampling Elbow, Filter	Adult	3 L Latex Free	20/ea
963-6754	UniLim, Gas Sampling Elbow/line	Pediatric	2 L	20/ea



Contact your Tri-anim Account Manager about other configurations to meet your anesthesia breathing circuit needs.

ICU Medical Anesthesia Breathing Circuits

ICU Medical's circuit systems deliver exceptional value by providing the configuration that meets your anesthesia department's specific needs. With several standard circuit configurations and a variety of components in either adult or pediatric sizes, regular or expandable tubing, as well as adult and pediatric single-limb circuits, we have a circuit that fits your requirements.

Item #	Description	Size	Bag	Qty
70-450954-NL	Dual limb, Expandable	Adult	3 L	20/ca
70-492804-NL	Expandable, 40-100 in, Elbow, Filter, Straight Wye	Pediatric	—	20/ca
70-650904	Expandable, 33-87 in, Elbow, Filter, Straight Wye	Adult	3 L	20/ca
70-370941	Corrugated, 40 in, Elbow, M/M Line, Wye, Filter, Mask	Adult	3 L	15/ca
70-372904	Corrugated, 60 in, Elbow, Straight Wye, Filter	Pediatric	1 L	15/ca
70-373600	Swivel Wye, Gas Port, 40 in, Filter	Adult	3 L	20/ca
70-370924	40 in, Filter	Adult	3 L	15/ca



70-370941

WE PARTNER WITH THE BEST, TO BRING YOU THE BEST.

Tri-anim partners with leading manufacturers to support anesthesiologists with new specialty anesthesia products and everyday consumables.



SPECIALTY ANESTHESIA PRODUCTS AND EVERYDAY CONSUMABLES

Tri-anim® offers products in anesthesia administration, airway management, monitoring and patient recovery. You can count on us to provide your anesthesia department with the right products when you need them.

THE TRI-ANIM ADVANTAGE



INNOVATIVE PRODUCTS



We partner with over 65 anesthesia product manufacturers



We offer more than 1,700 anesthesia products



CONSULTATIVE APPROACH



Over 100 clinically-trained Account Managers



Dedicated to your success from product identification to training and implementation



EXCEPTIONAL SERVICE



98% of all orders are delivered in 2 days or less



Over 40 in-house Customer Care representatives ready to answer your questions



Advancing Healthcare, Together.

Contact your Account Manager today!

▶ 800.874.2646

▶ TRI-ANIM.COM

customer account number