



# Curaplex® TritonGrip EC™

TritonGrip EC™ nitrile gloves provide excellent protection against abrasions, snags, and punctures, with chemical resistance to most liquids. These textured gloves feature an extended cuff and are 6-mil thick for extra protection. They are powder free, non-sterile and single use. No latex was used in the manufacturing of this product.

SPECIFICATIONS	
Material	Nitrile
Color	Blue
Length	12 in
Cuff thickness	6 mil
Palm thickness	5 mil
Finger thickness	6 mil

TENSILE STRENGTH	
Before aging	21 MPa min
After aging	21 MPa min

ELONGATION	
Before aging	500%
After aging	450%

## Ordering Information

- 121-47139      Small, 50/bx
- 121-47140      Medium, 50/bx
- 121-47141      Large, 50/bx
- 121-47142      X-Large, 50/bx

## Call for your quote!

For more information on pricing and promotions, please contact your Account Manager for your quote or call Customer Service at 800.558.6270

