

CURAPLEX® VINYL GLOVES



These gloves are versatile and convenient in a variety of settings - perfect for everything from food handling and preparation to janitorial applications. They are soft and pliable with a comfortable fit.

ORDERING INFORMATION

CUR5601	Small Vinyl Gloves	100/box
CUR5602	Medium Vinyl Gloves	100/box
CUR5603	Large Vinyl Gloves	100/box
CUR5604	X-Large Vinyl Gloves	100/box

SPECIFICATIONS

Material	Vinyl
Color	Clear
Length	9 in
Texture	Smooth
Palm thickness	.08 mm
Finger thickness	.08 mm
Hand Compatibility	Ambidextrous
Tensile Strength	13 MPa
Elongation at Break	300%
AQL	1.5

WEIGHT

S	4.5 g
M	5.0 g
L	5.5 g
XL	6.0 g

FEATURES

- Beaded Cuff
- Medical Grade
- Latex free
- Powder free
- 5 Year Shelf Life



800.558.6270

BUYEMP.COM



Emergency
Medical Products