

Pediatric Immobilization



Ferno's complete line of **pediatric immobilization** solutions provide **fast and safe spinal immobilization for pediatric patients** ranging in size from 2-90 lbs. These solutions are designed to work with adult-size backboards and ambulance cots so there's no need to purchase additional transport equipment.



Medkids® Pedi-Sleeve quickly adapts standard adult backboards for immobilization of patients from 12 lbs-60 lbs.



Pedi-Pac® Immobilizer provides fast, safe spinal immobilization and restraint for children from 28-54 inches tall, weighing 20-90 lbs.



Medkids® Baby Board safely secures infants ranging in size from 2-15 lbs.

Pediatric Immobilization

Medkids® Pedi-Sleeve



The Pedi-Sleeve quickly adapts standard adult backboards for pediatric immobilization of patients ranging in size from 12 lb (5 kg) to 60 lb (27 kg). Repositionable, crescent-shaped compressible head blocks, restraint straps, and a self-contained Pneumatic Positioning Device (PPD) ensure patients are secure and fully immobilized.

- Rolls into a self-contained carrying case to store in small compartments, drawers, or under bench seats
- Crescent-shaped compressible head blocks allow access and visual monitoring of the ears while maintaining immobilization
- Designed for infants and children ranging from 12 lb (5 kg) to 60 lb (27 kg)

Medkids® Baby Board



The Baby Board secures infants ranging in size from 2-15lb and features a self-contained Pneumatic Positioning Device (PPD) to help maintain cervical neutral alignment. The Baby Board is made of durable Dupont® Cordura® and is X-ray translucent.

- Compatible with most EMS cots and isolettes to provide safe and secure transport
- Crescent-shaped head blocks secure cervical spines for infants too small for standard size cervical collars
- Accommodates the wide range of thoracic elevation requirements without extra padding
- Unique design accommodates all central lines, umbilical lines, chest tubes, urinary catheters and other invasive and non-invasive equipment
- Designed for pre-term neonates and larger infants ranging from 2 lb (1 kg) to 15 lb (7 kg)

Pedi-Pac® Immobilizer



The Pedi-Pac® provides spinal immobilization and restraint for children from 28-54 in. tall and weighing 20-90 lb (9 to 41 kg). Cot-fastening loops, an adjustable head support, and color-coded straps provide fast, easy application and sewn-in lifting handles at both ends allow easy handling in confined areas.

- Adjustable head support with ear openings for monitoring fluid drainage
- Includes head and chin straps, and a carrying case
- Individual leg restraints allow one leg to be immobilized while a medic attends to the other leg

Learn more at www.fernoems.com/pediatricimmob