

NON-CONTACT INFRARED THERMOMETER

IR 200

This is an FDA 510k cleared device

Mode Indicator

**Easy To Understand
Digital Screen**

The screen background will illuminate to indicate no fever, slight fever, or high fever.

$\leq 99.2^{\circ}\text{F}$

No
Fever

$99.3^{\circ}\text{F} - 100.3^{\circ}\text{F}$

Slight
Fever

$\geq 100.4^{\circ}\text{F}$

High
Fever



Mode Button

Three temperature modes available: Body, Surface and Room.

Speaker Mode On/Off

Choose between audible mode or silent mode.

Memo Mode

View up to 32 past readings.



Infrared Sensor

Accurate for non-contact measurement.

On/Scan Button

Easy one-handed operation.

SPECIFICATIONS:

- **Measurement Time:** ≤ 1 second
- **Measuring Distance:** 1.2" - 2.0"
- **Body Temperature Range:** $89.6^{\circ}\text{F} - 109.4^{\circ}\text{F}$
- **Ambient Temperature:** $50^{\circ}\text{F} - 104^{\circ}\text{F}$
- **Accuracy:** $\pm 0.6^{\circ}\text{F}$
- **Automatic Power-Off:** < 30 seconds
- **Batteries:** 2 x AA
- **Product Longevity:** 100,000 uses
- **Product Weight:** 3.8 oz. (without batteries)
- **Product Size:** L 6.1" x W 3.6" x H 1.5"

ITEM #	DESCRIPTION	UOM
MS-131001	Non-Contact Infrared Thermometer, IR 200	1/EA, 40 EA/CS



800.558.6270

BUYEMP.COM



**Emergency
Medical Products**